

# Technical Data Sheet




**Product:**

**Tobaguard LS**

**Code:**

**CAT-TG10-LS**



## Description:

Tobaguard LS is a monitoring trap for *Lasioderma serricorne* beetles. The lure contains nature identical pheromones to attract the target insect to the trap where they will remain in the adhesive layer.

**Active Ingredient:** Serricorne pheromone.

**Inert Ingredient:** Craft board trap with adhesive coating.

**Longevity:** 4-6 weeks once opened.

**Packaging:** 10 Tobaguard traps with 10 Tobaguard rings.

## Application

Remove one trap from packaging. Fold the trap along the crease and peel off the protective paper. Open the foil sachet containing the lure. Do not open the ring capsule, apply capsule with cap facing down as shown in the image below. Assemble trap by inserting the tabs securely into the opposing slits. Place the trap on/near the wall in the desired location near to food commodity.

**Trap density:** One Tobaguard every 5-10 meters at a height of 1-1.5m.

**Trap data:** Examine the trap once a week, or at least once a month. During the height of the infestation, more frequent readings may be necessary.

## Storage

Store at room temperature in a cool, dry place.

**Longevity in storage:** 6-12 months, refrigeration is recommended to prolong shelf life.

**Replacement and disposal:** Replace every 6 weeks or as recommended by professional pest controller. Remove the used trap and dispose of away from the monitoring area once full with insects.

## Precautions

Wash hands before and after contact.

If irritation occurs, seek medical advice.